

## PRODUCT SPECIFICATION CHANGE

July, 2009

Dear Valued Customer:

Data compiled from a review of product lot results for the following product has led to a change in the osmolality specification:

**Description:** RPMI 1640 with L-glutamine & 25mM HEPES  
**Catalog Number:** 10-041-CV, 10-041-CM

Beginning with lot 10041098, scheduled to be produced in September 2009, the osmolality specification will be changed from 260 – 320 mOsm/kg H<sub>2</sub>O to 275 – 315 mOsm/kg H<sub>2</sub>O. There have been no changes made to the formulation. This change is to reflect a more accurate representation of the osmolality value for this product that Mediatech has achieved over the years.

Mediatech makes every effort to ensure the integrity and quality of our products before and after they leave our facility. The on-site product specification review program enables Mediatech to monitor test results and overall product integrity on a continuous basis and make adjustments as needed.

If you have any questions or concerns regarding this specification change, please do not hesitate to contact our Technical Support department.

1.800.CELLGRO ..... TOLL FREE  
703.471.5955 ..... LOCAL AREA  
703.471.1086 ..... FACSIMILE