

Change Notification July, 2010

Dear Valued Customer:

The purpose of this letter is to inform customers of a Raw Material Change to film used in the manufacture of some flexible containers and tank liners offered by Mediatech, Inc.

Mediatech has been notified by our supplier of a discontinuation/replacement of an anti-blocking agent used in the manufacture of a Polyethylene film to make some bags and tank liners. The new raw material has been compared technically to the raw material being replaced, and has been found to be a comparable and acceptable replacement from a technical standpoint. The new material is animal free. Further more the Drug Master File (DMF # 8686), Certifications and Final Specifications for this film, will not change.

The supplier anticipates this change to take place for orders occurring in and after late July to early August 2010.

The following Product numbers will be impacted by this change:

PART	DESCRIPTION
91-200-50	50 L BAG
91-200-62	100 L GUSSETED BAG
91-200-72	100 L BAG

1.800.CELLGRO TOLL FREE
703.471.5955 LOCAL AREA
703.471.1086 FACSIMILE

PART	DESCRIPTION
91-300-15	50 L GUSSETED TANK LINER
91-300-20	130 L TANK LINER
91-300-25	100 L TANK LINER
91-300-30	200 L TANK LINER
91-300-50	400 L TANK LINER
91-300-60	560 L TANK LINER
91-300-70	790 L TANK LINER
91-300-80	1090 L TANK LINER
91-300-90	1360 L TANK LINER
91-600-67	2 L BAG
91-600-71	20 L BAG W/POUR SPOUT

Mediatech makes every effort to ensure the integrity and quality of our products before and after they leave our facility. The on-site product specification review program enables Mediatech to monitor test results and overall product integrity on a continuous basis and make adjustments as needed.

If you have any questions or concerns regarding this specification change, please do not hesitate to contact our Technical Support department.

1.800.CELLGRO	TOLL FREE
703.471.5955	LOCAL AREA
703.471.1086	FACSIMILE