

Change Notification

April, 2010

Dear Valued Customer:

Effective January 1, 2010 there was a change in a polypropylene component used in the manufacture of Clear and Opaque C-Flex® formulations. C-Flex® tubing is manufactured by Saint-Gobain Performance Plastics/ Clearwater and supplied to Mediatech, Inc as components to bag assemblies. This change was due to a supplier discontinuation of the polypropylene component. Saint-Gobain has selected a replacement polypropylene and validated it for ASTM physical property testing, biocompatibility testing, extruding and moldability, sealing and welding capabilities, pump life and shelf life testing. Validation testing was conducted on both pre-sterilized and post-sterilized samples of C-Flex® resins. As a result of the polypropylene change, all C-Flex® Clear and Opaque formulations are now Animal Derived Component Free.

The Master Validation Plans are available to assist with answering any questions or concerns you may have.

The following Standard Catalog Items are affected by this change:

Part No.	Description
91-100-30	50L Collection Container
91-100-35	100L Collection Container
91-200-01	1L Collection Container, Gamma Irradiated
91-200-02	2L Collection Container, Gamma Irradiated
91-200-05	5 L Collection Container, Gamma Irradiated
91-200-10	10 L Collection Container, Gamma Irradiated
91-200-20	20L Collection Container, Gamma Irradiated
91-200-50	50 Liter Gusseted Bag
91-700-12	Filling Bell Assembly

Mediatech makes every effort to ensure the integrity and quality of our products before and after they leave our facility. The on-site product specification review program enables Mediatech to monitor test results and overall product integrity on a continuous basis and make adjustments as needed.

If you have any questions or concerns regarding this specification change, please do not hesitate to contact our Technical Support department.

1.800.CELLGRO TOLL FREE
 703.471.5955 LOCAL AREA
 703.471.1086 FACSIMILE