

Change Notification

April, 2010

Dear Valued Customer:

Mediatech is making a raw material change to the following products:

Catalog Number	Description	Present formulation with EDTA tetrasodium salt tetrahydrate	Formulation change with EDTA tetrasodium salt dihydrate
25-051-CI	0.05% Trypsin / 0.53 mM EDTA in HBSS without calcium and magnesium	0.24 g/L	0.221 g/L
25-052-CI 25-052-CV	0.05% Trypsin / 0.53 mM EDTA in HBSS without sodium bicarbonate, calcium and magnesium	0.24 g/L	0.221 g/L
25-053-CI	0.25% Trypsin / 2.21 mM EDTA in HBSS without sodium bicarbonate, calcium and magnesium	1.0 g/L	0.920 g/L

Mediatech can no longer source EDTA, tetrasodium salt, tetrahydrate in a consistent manner from approved suppliers so the above formulations will now use EDTA, tetrasodium salt, dihydrate with adjustments as indicated in the above table. In addition, two proprietary formulations are affected as well. All formulations are being changed so as to maintain the molar concentrations of EDTA and will not have any impact on the quality of the solutions.

Mediatech makes every effort to ensure the integrity and quality of our products before and after they leave our facility. The on-site product specification review program enables Mediatech to monitor test results and overall product integrity on a continuous basis and make adjustments as needed.

If you have any questions or concerns regarding this specification change, please do not hesitate to contact our Technical Support department.

1.800.CELLGRO TOLL FREE
 703.471.5955 LOCAL AREA
 703.471.1086 FACSIMILE