

Change Notification

May, 2010

Dear Valued Customer:

Mediatech is making a raw material change to the following product:

Catalog Number	Description	Present formulation with yeastolate	Formulation change with ultra-filtered yeast extract
13-100-CV	Hink's TNM-FH Medium (Mod.), 1X	3.333 g/L	3.333 g/L

Ultra-filtered yeast extract is a more pure grade of material and formulations using it performed as well or better than formulations with yeastolate. The solution will continue to have the same specifications for physiochemical properties.

Mediatech makes every effort to ensure the integrity and quality of our products before and after they leave our facility. The on-site product specification review program enables Mediatech to monitor test results and overall product integrity on a continuous basis and make adjustments as needed.

If you have any questions or concerns regarding this specification change, please do not hesitate to contact our Technical Support department.

1.800.CELLGRO TOLL FREE
703.471.5955 LOCAL AREA
703.471.1086 FACSIMILE