



# HAT (Hypoxanthine, Aminopterin, Thymidine) HT (Hypoxanthine, Thymidine) Solutions

Catalog No.	25-046 50x, Liquid HAT	25-047 50x, Liquid HT
Components	mg/L	mg/L
Aminopterin	8.81	--
Hypoxanthine, Na	910.70	910.70
Thymidine	193.80	193.80
Specifications		
pH (before adjustment)	--	--
Final pH	10.5 ± 0.5	--
Osmolality (mOsm)	9 ± 15	12 ± 10