



## MEM Amino Acids Solution

Catalog No.	25-030 50X Liquid MEM Amino Acid Solution
Components	mg/L
L-Arginine • HCl	6320.00
L-Cystine • 2HCl	1560.00
L-Histidine • HCl • H <sub>2</sub> O	2095.00
L-Isoleucine	2625.00
L-Leucine	2625.00
L-Lysine • HCl	3625.00
L-Methionine	750.00
L-Phenylalanine	1625.00
L-Threonine	2380.00
L-Tryptophan	500.00
L-Tyrosine • HCl	2160.00
L-Valine	2340.00
Specifications	
pH	1.5 ± 0.5
Osmolality (mOsm/kg)	475 ± 23