



MEM Nonessential Amino Acids Solution (NEAA)

Catalog No.	25-025 100X Liquid MEM Nonessential Amino Acid Solution
Components	mg/L
L-Alanine	890.00
L-Asparagine • H ₂ O	1500.00
L-Aspartic acid	1330.00
L-Glutamic acid	1470.00
Glycine	750.00
L-Proline	1150.00
L-Serine	1050.00
Specifications	
pH	3.5 ± 0.5
Osmolality (mOsm/kg)	74 ± 15