



MEM Vitamin Solution

Catalog No.	25-020 100X Liquid MEM Vitamin Solution
Components	mg/L
NaCl	8500.00
D-Calcium pantothenate	100.00
Choline Chloride	100.00
Folic acid	100.00
i-Inositol	200.00
Nicotinamide	100.00
Pyridoxine • HCl	100.00
Riboflavin	10.00
Thiamine • HCl	100.00
Specifications	
pH	5.8 ± 1.3
Osmolality (mOsm/kg)	250 ± 50