



Penicillin-Streptomycin Solutions

Catalog No.	30-001 Penicillin-Streptomycin 50X	30-002 Penicillin-Streptomycin 100X	30-009 Penicillin-Streptomycin- L-Glutamine
Components	g/L	g/L	g/L
K Penicillin G	5,000,000 U/L	10,000,000 U/L	10,000,000 U/L
Dihydro-Streptomycin • SO ₄	5.0	10.0	10.0
L-Glutamine	--	--	29.2
Na Citrate • 2H ₂ O	--	--	2.582
NaCl	--	--	1.4
Specifications			
pH	INFO	INFO	6.0 ± 0.5
Osmolality (mOsm/kg)	INFO	INFO	350 ± 50