

Trace Elements Solutions

Catalog No.	99-182 1000X Liquid Trace Elements A
Components	mg/L
CuSO ₄ • 5H ₂ O	1.60
ZnSO ₄ • 7H ₂ O	863.00
Selenite • 2Na	17.30
Ferric citrate	1155.10
Specifications	
pH	6.25 ± 0.25
Osmolality (mOsm/kg)	INFO

Catalog No.	99-175 1000X Liquid Trace Elements B
Components	mg/L
MnSO ₄ • H ₂ O	0.17
Na ₂ SiO ₃ • 9H ₂ O	140.00
Molybdic acid, Ammonium salt	1.24
NH ₄ VO ₃	0.65
NiSO ₄ • 6H ₂ O	0.13
SnCl ₂ (anhydrous)	0.12
Specifications	
pH	INFO
Osmolality (mOsm/kg)	INFO

Catalog No.	99-176 1000X Liquid Trace Elements C
Components	mg/L
AlCl ₃ • 6H ₂ O	1.20
AgNO ₃	0.17
Ba(C ₂ H ₃ O ₂) ₂	2.55
KBr	0.12
CdCl ₂	2.28
CoCl ₂ • 6H ₂ O	2.38
CrCl ₃ (anhydrous)	0.32
NaF	4.20
GeO ₂	0.53
KI	0.17
RbCl	1.21
ZrOCl ₂ • 8H ₂ O	3.22
Specifications	
pH	INFO
Osmolality (mOsm/kg)	INFO