

## L-Glutamine Addition Chart for cellgro® Media

Media	mg/L	mM	mL/L*
AMEM	292.00	2.0	10
BME	292.30	2.0	10
DMEM	584.00	4.0	20
DMEM/F12	362.10	2.5	15
F12K	292.00	2.0	10
Ham's F-10	146.2	1.0	5
Ham's F-12	146.2	1.0	5
Iscove's DMEM	584.00	4.0	20
M199	100.00	0.68	1.7
EMEM	292.00	2.0	10
RPMI 1640	300.00	2.05	10

Product Description	Catalog No.	Size
L-Glutamine Powder, MW 146.15	61-030-RM	1 x 100 g
L-Glutamine 200 mM Solution	25-005-CI 25-055-CV	6 x 100 mL 6 x 500 mL