

Sodium Bicarbonate Addition Chart for cellgro[®] Powder Media

Media	Catalog No.	g/L	mL/L*
AMEM	50-012	2.2	29.3
	90-012	2.2	29.3
BME	50-005	2.2	29.3
	50-006	2.2	29.3
DMEM	50-003	3.7	49.3
	50-013	3.7	49.3
	50-014	3.7	49.3
	90-013	3.7	49.3
	90-113	3.7	49.3
DMEM/F12	90-090	2.438	32.0
	90-091	2.438	32.0
Ham's F10	50-030	1.2	16.00
Ham's F12	50-040	1.2	16.00
Iscove's DMEM	50-016	3.0	40.0
	90-016	3.0	40.0
Leibovitz's L-15	50-045	0.0	0.000
M199	50-050	2.2	29.3
	50-051	2.2	29.3
	90-050	2.2	29.3
EMEM	50-004	2.7	37.3
	50-010	2.2	29.3
	50-011	2.2	29.3
	50-019	2.2	29.3
	51-010	2.2	29.3
	90-009	2.2	29.3
McCoy's 5A	50-025	2.2	29.3
RPMI 1640	50-020	2.0	26.7
	90-022	2.0	26.7

Product Description	Catalog No.	Size
Sodium Bicarbonate Powder, MW 84.01	61-065-RO	1 x 500 g
* 7.5% Solution	25-035-CI	6 x 100 mL